



SAFESEPT^{MAX} Handdisinfection Gel

- alcohol based



PRODUCT DESCRIPTION:

SAFESEPT MAX[®] Handdisinfection Gel is free from cumulative long-term acting ingredients that might have an irritating or sensitising effect on the skin. Precise dosing and easier, non-drip application through a user friendly gel formulation. Compatible with all washing lotions. Free of aldehydes, phenols and quaternary ammonium compounds.

AREAS OF APPLICATION:

Ready-to-use alcohol disinfectant for hygienic hand disinfection.

SPECTRUM OF ACTION:

Bactericidal, yeasticidal, tuberculocidal, virucidal (HBV, HIV, HCV).

APPLICATION TIME:

Hygienic hand disinfection (EN 1500): 2x (3ml – 30 secs.)

APPLICATION:

For hygienic hand disinfection, rub hands with undiluted solution and keep them damp for 30 seconds

COMPOSITION:

100 g contains: 63g ethanol (CAS: 64-17-5)

SPECIAL INFORMATION:

Highly flammable liquid and vapour. Causes serious eye irritation. Harmful to aquatic life with long lasting effects. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Dispose of contents/container to approved disposal company or local collection. For Professional Use Only.





SAFESEPT^{MAX} Handdisinfection Gel

Handdisinfection Gel

- alcohol based



PRODUCT STATUS/LISTINGS:

- Conform to the EU Biocidal Product Regulation (BPR, 528/2012).

AVAILABLE DELIVERY FORMS:

- | | |
|--------------|----------------------------|
| REF 988-2737 | 300ml dispenser bottle |
| REF 988-2738 | 500ml Eurodispenser Bottle |
| REF 988-2756 | 1L Eurodispenser Bottle |
| REF 988-2757 | 5L Canister |

AVAILABLE ACCESSORIES:

- | | |
|--------------|-----------------------------------|
| REF 979-6274 | Dosing pump for 500ml |
| REF 979-6275 | Dosing pump for 1L |
| REF 979-6276 | Tab/dosing device for 5L canister |
| REF 979-6278 | Canister-key |
| REF 979-6272 | Dosing pump for 5L Canister |